

## ChIP-seq protocol

- Grow 2ml overnight culture.
- Next morning culture- subculture 1/100 (e.g. 400µl into 40ml). Grow to appropriate OD600 (0.3-0.6).
- Cross-link cells by adding formaldehyde to a final concentration of 1% (1ml of 37% stock to a 40ml culture). Swirl to mix and leave for 20 minutes.
- Quench formaldehyde by adding glycine to a final concentration to 0.5M (typically 10ml of 2.5M glycine per 40ml culture).

**NOTE: if using a smaller culture than 40ml, change the amount of formaldehyde used. Keep the volume of glycine the same.**

- Spin cells down for about 5 minutes at 4000rpm. Re-suspend in TBS (0.5-1 times the culture volume- resuspend in 1ml first, then add the rest)). Spin down again.
- Resuspend in 1ml TBS and transfer to eppendorf. Spin in microcentrifuge at max speed for minute.
- Discard supernatant. Cells can be frozen at this point.
- Resuspend cells in 1ml FA lysis buffer (150mM NaCl) containing 4mg/ml lysosome. Incubate at 37°C for 30 minutes to lyse cells. Crosslinked lysates can be frozen at this point.
- Chill on ice for >5 minutes. Sonicate as directed:

**USING A TIP SONICATOR: 2 x 30s pulses, 50% output, keeping cells on ice.**

**USING A BIO-RUPTOR: Turn on pumps at least 45 minutes before use to allow water bath to chill to 4°C. Optimise pump speed to maintain marked water level. 2 x 15 minute cycles (check water level in between).**

- Spin cell lysates in microcentrifuge at max speed for 5 minutes. Keep supernatant (about 1ml). If cultures >20ml in volume dilute in FA lysis buffer (150mM NaCl) so there is 1ml of lysate for each 20ml of original culture. (E.g. 40ml culture- dilute to 2ml).
- Set up cocktails containing 500µl lysate and 300µl FA lysis buffer (150mM NaCl). (4 per 40ml culture). Use 1 of these for each immunoprecipitation.
- OPTIONAL: save 20µl of lysate as 'input'. De-crosslink before you do anything with it.
- Cross-linked and sonicated lysates can be frozen at this point.
- Immunoprecipitation:

#### IMMUNOPRECIPITATION:

- Use one of the 800µl cocktails as described in (11). i.e. 10ml culture sonicated lysate.
  - Add 25µl of Protein A beads (50% slurry with TBS). Resuspend gently before use. BLUNT TIPS before use (use scissors).
  - Add 1µl of Neoclone<sup>®</sup> antibody β/σ/CRP or 2µl of anti-FLAG (Sigma<sup>®</sup>) anti-H-NS.
  - Rotate at room temperature for 90 minutes. Parafilm lid. For CRP, rotate at 4°C overnight.
- Spin at 4000rpm for 1 minute, remove the supernatant carefully. Add 700µl of FA lysis buffer (150mM NaCl) and resuspend gently using a blunt tip. Transfer to a spin-X-column.
  - Rotate for 3 minutes at room temp (check lids are closed before rotating), then spin at 4000rpm for 1 minute. Discard supernatant.

WASHES IN SPIN-X-COLUMN: With all wash steps, use 750µl unless indicated; check lids are closed and rotate at room temp for 3 minutes before spinning at 4000rpm for 1 minute. Discard supernatant unless directed otherwise. It is useful after a set of washes to briefly re-spin the spin-X-columns for a minute to remove any excess liquid left.

- Washes:
  - FA lysis buffer 150mM NaCl.
  - FA lysis buffer 150mM NaCl.
  - 10mM Tris-HCl, pH7.5.
  - 10mM Tris-HCl, pH7.5.

- Blunting: use NEB Blunting kit:

**BLUNTING USING NEB® KIT (resuspend in buffer then add enzyme). Add:**

- 10µl 10x quick blunting buffer.
- 10µl dNTP mix (supplied in kit).
- 80µl dH<sub>2</sub>O.
- 2µl blunt enzyme kit (T4 DNA polymerase).

Parafilm lid and rotate (in such a way that doesn't invert) for 30 minutes at room temp.

- Washes:
  - FA lysis buffer 150mM NaCl.
  - FA lysis buffer 150mM NaCl.
  - 10mM Tris-HCl, pH8.
  - 10mM Tris-HCl, pH8.
- A tail: use NEB Klenow (3' -> 5' exo-) kit:

**Resuspend beads in the following mixture:**

- 10µl 10x NEB buffer 2
- 2µl 100mM dATP (supplied).
- 88 µl dH<sub>2</sub>O.
- 2µl Klenow.

Parafilm lid. Incubate at 37°C rotating (not inverting) for 30mins).

- Washes:
  - FA lysis buffer 150mM NaCl.
  - FA lysis buffer 150mM NaCl.
  - 10mM Tris-HCl, pH7.5.
  - 10mM Tris-HCl, pH7.5.

- Ligate adaptors:

**Resuspend beads in the following mixture:**

- 100µl 1x ligase buffer (50µl dH<sub>2</sub>O + 50µl 2x ligase buffer).
- 1µl NEXTflex ChIP-seq barcoded adaptor (BioO) (have used 1µl non-diluted before).
- 4µl quick ligase.

**Parafilm lid and rotate (not invert) at room temp for 15 mins.**

- Washes:
  - FA lysis buffer 150mM NaCl.
  - FA lysis buffer 150mM NaCl.
  - FA lysis buffer 500mM NaCl.
  - ChIP wash buffer.
  - TE
  - Flush 300µl TE through the column after just to get rid of any bubbling.
- Transfer spin-X-columns to new dolphin-nosed tubes, add 100µl elution buffer and incubate at 65°C for 10 minutes.
- Elute by quickly transferring to a microcentrifuge: spin at 4000rpm for 1 minute, transfer flowthrough to an eppendorf.
- To de-crosslink: boil for 10 minutes, use a tube lock if possible. De-crosslink inputs at this point as well, if used.

- Phenol extract OR use Ampure® magnetic beads to clean up.

#### To phenol extract:

- Add 300 $\mu$ l dH<sub>2</sub>O and 400 $\mu$ l PCIA (pH6.7/8.0), vortex.
- Take upper phase and transfer to new tube.
- Add 1 $\mu$ l 35 mg/ml glycogen, 40 $\mu$ l 3M Na acetate, and 1ml ethanol.
- Vortex and place at -20°C for > 2hours.
- Spin at max speed at 4°C for 15-20mins. Longer is better.
- Dump supernatant, add 1ml 80% ethanol.
- Spin again for 15mins at max speed at 4°C.
- Dump supernatant. Spin briefly again and remove any residual ethanol.
- Dry, either 10mins at 65°C with lids open or speedvac.

#### Ampure® magnetic beads:

- Warm up beads to room temperature, and resuspend gently.
- Add 1.1x volume beads to DNA solution.
- Mix gently 10x (blunt tips helpful), leave at room temperature for 5 minutes.
- Place on magnet for 2 minutes.
- Remove solution from beads.
- KEEP ON MAGNET. Add 200 $\mu$ l 70% ethanol for 30 seconds, and then remove ethanol. REPEAT.
- Air-dry for at least 5 minutes (65°C for 10 minutes works better).
- Off the magnet, resuspend beads in 11/12 $\mu$ l of water.
- Place on magnet for 1 minute, remove and keep liquid.

#### Amplifying library

- First run a real time- PCR, use 1 $\mu$ l of ChIP DNA as template- dilute 1/8, use 2 $\mu$ l each in three wells.
- OLIGOS= JW3764/3765 (primer mix for NEXTflex, 10 $\mu$ M each, 3 $\mu$ l of this mix per reaction.)
- Analyse on manual Ct set to 0.1.
- Cycle number at which sample 'appears'/ 'comes up'= X.

- Add 3 to X and run a PCR with this number of cycles. If using multiple samples, open up PCR machine after elongation cycle to take out tubes.
- Each mix (set up 4 reactions in strip PCR tubes):
  - 5µl 10x taq buffer.
  - 2µl DNA with ligated adaptors.
  - 1µl JW3764 100µM.
  - 1µl JW3765 100µM.
  - 0.4µl dNTP mix (25mM each).
  - 1.5µl taq.
  - 40.1µl dH<sub>2</sub>O.
  - TOTAL= 50µL.

#### Clean and concentrate library

- Clean up with Ampure® magnetic beads, QIAGEN MinElute®, or ZymoGen clean and concentrate®. Elute in 15-25µl dH<sub>2</sub>O (15µl better- then add 10µl loading dye.)
- Load onto a non-denaturing 8% acrylamide gel.

#### **8% acrylamide gel:**

- 7ml dH<sub>2</sub>O.
- 1ml 10xTBE.
- 2ml 40% acrylamide.
- 100µl 10% APS.
- 10µl TEMED.

- Run for 1 hour at 100V.
- Stain in Ethidium bromide for 5 minutes.
- On trans-illuminator, cut out a smear from 200-600bp, avoid the 150bp band- this is the adaptor-adaptor band.
- Pierce the bottom of a 0.7ml tube with a hot needle, place inside a 1.7ml tube.
- Place gel slice into 0.7ml tube and spin at max speed for 5 minutes to shred the gel.
- Remove 0.7ml tube, add 0.4ml acrylamide gel extraction buffer, parafilm lid and rotate overnight.
- Transfer gel and liquid to a spin-X-column (blunt tips helpful for removing the gel).
- Spin at max speed for 5 minutes.
- Precipitate DNA as follows:

**Precipitation:**

**Add:**

- 1ul 35mg/ml glycogen.
- 0.4ml propan-2-ol/ 1ml cold ethanol.

**Vortex and place at -20°C for at least 2 hours.**

- Spin for 15-30 mins at max speed at 4°C, dump supernatant.
- Add 1ml 80% ethanol, spin for 10-15 mins again.
- Remove supernatant, spin for 3-4mins, remove any remaining supernatant.
- Dry at 65°C for 10 minutes.
- Resuspend in 11-12µl dH<sub>2</sub>O.

**Quantify library:**

- Use a Qubit® 2.0 Fluorometer - make up working solution in the following ratio:
  - 1µl 200x dye : 199µl dilution buffer.
- Dilute standards 1 & 2 as follows:
  - 10µl standard : 190µl working solution.
- Dilute samples 1µl in 199µl of working solution.
- MAKE UP ALL TUBES IN THE THIN-WALLED TUBES SPECIFIED.
- 0.5µg/L to 50µg/L is expected (when diluted 1:200).
- Multiply by 200, divide by 1000 to get ng/µl.
- Collate samples in an equimolar ratio, using the whole volume of the most dilute solution. Can concentrate sample in speedvac, 10µl works well (sent 6µl of this off for sequencing).

**CAT. NUMBERS**

Monoclonal Anti-FLAG (R) M2, antibody PR&	SIGMA	F3165-1MG
Spin-X .2CA, 25/BAG-N/S	VWR	83-8161-CS
2.0ml dolphin nosed tubes 250/bag Corning	VWR	29442-590
Quick Blunting™ and Quick Ligation™ Kits	NEB	E0542L
Klenow Fragment (3'→5' exo-)	NEB	M0212L
Phase Lock tubes GEL HVY 2ML 200PK	Fisher Sci	FP2302830
NEXTflex ChIP-seq Barcodes- 24	Bio0 Scientific	514122
Agencourt AMPure XP 5 ml	Beckman Coulter	A63880

kit, 139/278 preps		
DynaMag -2Magnet	Life Technologies	12321D
Qubit Fluorometer 2.0	Life Technologies	Q32866
Qubit® dsDNA HS assay Kit (for use with Fluorometer)	Life Technologies	Q32851

**Illumina primer sequences (for library amplification)**

AATGATACGGCGACCACCGAGATCTACAC

NEXTflex CHIP Primer 1

CAAGCAGAAGACGGCATACGAGAT

NEXTflex CHIP Primer 2



## Buffers

### FA lysis buffer

50mM Hepes-KOH, pH7

150/500mM NaCl

1mM EDTA

1% Triton X-100

0.1% Sodium deoxycholate

0.1% SDS

### ChIP wash buffer

10mM Tris-HCl, pH 8.0

250mM LiCl

1mM EDTA

0.5% Nonidet-P40

0.5% Sodium Deoxycholate

### ChIP elution buffer

50mM Tris-HCl, pH 7.5

10mM EDTA

1% SDS

### TE

10mM Tris-HCl, pH 8.0

1mM EDTA

### 1x TBS

20mM Tris-HCl, pH7.4

0.9% NaCl