

FIS purification

tartar grew cultures in Amp+Kan but strain he gave me is Amp+Tet?! => Had given me wrong strain, replaced with correct one.

- spermidine 85559-1G £18.80 sigma
- lysozyme C7651 10g £134.00 sigma
- Heparin sepharose 17-0998-01 £270.00 GE Healthcare
- Heparin sepharose columns 5x5ml 17-0407-03 £379.00 GE Healthcare
- Resource S column 17-1180-01 £1,050.00 GE Healthcare

↳ (in the end used a High trap SP-XL cation exchange column)

- FIS test overexpression. 22/7/08
- Inoculate ^{2x} 10ml of LB (+amp+kan) with 200µl of O/N RJ4529 culture.
- When cultures get to an OD₆₅₀ of ~0.6 add 12.5µl of 1M IPTG and grow ^{both} for a further 60 minutes.
to one culture

14/7/08

- Use 15ml of O/N CB culture of RJ4529 (with Amp+Kan) to inoculate 750ml of pre-warmed LB.
- Grow to and OD₆₅₀ of ~0.6
- Add 940µl of 1M IPTG. Continue growing @ 37°C for a further 1 hour.
- Spin down cells and resuspend in 15ml 50mM Tris-HCl (pH8) + 10% Sucrose
- Freeze cells in -80°C then thaw.
- Add:
 - 71µl 0.5M EDTA
 - 19µl 0.1M PMSF
 - 94µl 1M DTT
 - 94µl 1g/ml spermidine
 - 1.9µl 3M Ammonium SO₄²⁻
 - 38µl 100mg/ml lysozyme

⇒ Incubate on ice for 90 minutes, then sonicate ^{3x30 sec}

⇒ Spin down debris at high speed and keep supernatant (can freeze)

FB buffer: 4L of ... 20mM Tris-HCl pH7.5
 of 20mM Tris-HCl pH7.5
 1.0mM EDTA
 10% Glycerol
 1.5M NaCl

1.0mM EDTA
 10% Glycerol
 0.3M NaCl

15/7/08

- Thaw supernatant from yesterday and transfer to dialysis membrane.
- ~~Dialyze against 1L FB @ 4°C for 3hrs.~~
- ~~Change to a fresh 1L FB buffer and dialyze for further 90min.~~
 ⇒ Didn't do this, dialyzed against 2L FB O/N instead.


Column ①: 5ml Heparin Sepharose equilibrate with buffer FB.

16/7/08

- Load crude extract.
- Elute in gradient upto 1.5M NaCl. Should come off ~0.75M NaCl.

Gel (loaded 6µl blue + 4µl protein)

- 1 Marker
- 2 Crude lysate
- 3 ~~Flow through (F3)~~ Dialyzed lysate
- 4 ~~F23~~ F3
- 5 F25
- 6 F27
- 7 F29
- 8 F31
- 9 F33
- 10 F35
- 11 F37
- 12 F39
- 13 F41
- 14 F43
- 15 F45



⇒ FRACTIONS 40 → 46 STORED @ -80°C

ification re-run

Overexpress etc as previously done.

- After o/p dialysis centrifuge to remove any precipitate.

Prepare Akta

- Remove any column if attached.
- Set to load and run 100% B through line B.
- Run 100% A (also on load) through line A.
- Set to inject, switch lines on superloop, and wash out superloop.
- Switch loop lines again and wash out in opposite direction.
- Turn down flow rate, ^(2 ml/min) attach column, and wash to with A.
- While still on load, add sample to superloop. (pause?)
- Switch to inject and start collecting 4ml fractions.
- When protein is loaded, switch back to load.
- Start collecting 2ml samples.
- Set gradient to run over 180ml upto 100% B.
- Pool peak fractions 60 → 64 (10ml total in ~ 864ml NaCl)
- Make 1L buffer FB with 0.2M NaCl to use as buffer A with cation exchange column. ^(1ml Hi Trap SP-XL column) (20mM Tris, 1.0mM EDTA, 10% glycerol, 0.2M NaCl)
- Make 50ml FB with no salt to dilute peak fractions. → dilute 1:4 to give a 200mM NaCl


iel

- 1 M
- 2 Crude extract
- 3 Dialysed
- 4 Pooled Hep-Seph frac^{ns}

- 5 F 24
- 6 F 25
- 7 F 26
- 8 26 F 27
- 9 27 F 28
- 10 28 F 29
- 11 29 F 30
- 12 30 F 31
- 13 31 F 32
- 14 32 F 33
- 15 33 F 34

From H. Trap SP-XL (1ml) column

- Each lane
4µl protein +
6µl lye. loaded
5µl of this.



F34
F33
F32
F31
F30
F29
F28
F27
F26
F25
F24

Keep these fractions

Pooled + diluted frac. vs
from Hep-sepharose column
Dialyzed lysate
Crude lysate

M

- Used 5µl for each frac. in a Biorad assay.

Fraction	OD ₅₄₅	[Protein]
F28	0.040	10.2 µM
F29	0.065	16.6 µM*
F30	0.052	13.3 µM
F31	0.030	7.7 µM

⇒ Molecular weight of Fis is 11.24 kD.

* Tested in EMSA @ djs vs cbpA. Behaves as previous prep from Martin