

Primer extension

DEPC → use at 0.1%  
v/v

(A) Isolation of RNA

- ✓ 1. From a fresh transformation plate or restreak of your strain, grow up an overnight culture in 5 ml LB (plus antibiotic).
- ✓ 2. Inoculate 10 ml LB (plus antibiotic) with 200  $\mu$ l overnight culture, and grow culture to OD<sub>650</sub> of 0.4-0.5 (approx 2-3 h).
- ✓ 3. ~~Mix 2.5 ml culture with 4.5 ml RNAlater (Ambion) in 15 ml centrifuge tube. Prepare three of these tubes per culture.~~
- ✓ 4. ~~Incubate at room temperature for at least 5 minutes,~~ then spin for 20 min, 4000 rpm, 20°C. Remove supernatant. Cell pellets can be used immediately for the RNA preps, or can be stored at -20/-70°C.
- ✓ 5. Resuspend all three pellets in a total of 200  $\mu$ l TE buffer + 40  $\mu$ g/ml lysozyme. (i.e. resuspend first pellet in 200  $\mu$ l TE, then use this suspension to resuspend second pellet etc.)  
*FRESH per reaction (1  $\mu$ l/ml)*
- ✓ 6. Incubate for 15 min at room temperature.
- ✓ 7. Prepare buffer RLT by adding 10  $\mu$ l  $\beta$ -mercaptoethanol per 1 ml buffer RLT. Add 700  $\mu$ l buffer RLT, vortex, then add 500  $\mu$ l ethanol to each sample and mix by pipetting.
- ✓ 8. Split each sample between two RNeasy mini columns. Centrifuge 15 s at 10,000 rpm.
- ✓ 9. Discard flow-through and add 700  $\mu$ l buffer RW1 to column. Centrifuge 15 s at 10,000 rpm. Discard flow-through and place column into new 2 ml collection tube.
- ✓ 10. Add 500  $\mu$ l buffer RPE to column. Centrifuge 15 s at 10,000 rpm. *NU tube*
- ✓ 11. Add a further 500  $\mu$ l buffer RPE to column and centrifuge 2 min at 10,000 rpm.
- ✓ 12. Place column into new collection tube. Centrifuge at full speed for 1 min to remove residual ethanol.
- ✓ 13. Transfer column to new microfuge tube. Elute RNA from each column in 30  $\mu$ l RNase-free water. Combine the two eluates for each sample. You should have approx 60  $\mu$ l in total for each RNA prep. Keep RNA samples on ice at all times. *let stand for 20s*
- ✓ 14. To remove DNA contamination:
  - ✓ a. To each RNA sample, add:
    - 3.5  $\mu$ l 10x Turbo DNase Buffer
    - 1  $\mu$ l Turbo DNase (Ambion)
  - ✓ b. Incubate 30 min at 37°C.
  - ✓ c. Add 6.1  $\mu$ l DNase inactivating reagent and incubate at room temperature for 2 min, occasionally tipping/flicking the tube to mix.

d. Centrifuge 1.5 min at 10,000 rpm. Transfer supernatant to a clean tube and keep on ice.

- ✓ 15. Calculate RNA concentration using quantifier, and run on gel to check integrity of RNA (you should see two clear bands, corresponding to 16 s and 23 s rRNA). Make up the agarose and 1 x TBE in DEPC/treated water, and wash the tank thoroughly with Ethanol before use.
- ✓ 16. Freeze RNA in aliquots. It is advisable to aliquot the RNA such that each tube contains the required amount of RNA for a single primer extension reaction (20-30 µg, see section (C) below).

### (B) Labelling of primer

- ✓ 1. Mix the following in an RNase-free microfuge tube:

0.5 µl 100 µM primer (see Note 1)

2 µl 10 x polynucleotide kinase buffer

15.5 µl DEPC-treated water

1 µl  $\gamma^{32}\text{P}$ -ATP

1 µl T4 polynucleotide kinase (NEB)

20 µl

- ✓ 2. Incubate 30 min at 37°C, then 10 min at 68°C to inactivate the enzyme.
- ✓ 3. Store labelled primer at -20°C.

### (C) Primer annealing/extension

#### Day 1

- ✓ x 1. Precool centrifuge to 4°C and set heating blocks to 50°C and 75°C.
- ✓ 2. Mix 1 µl labelled primer with 20-30 µg RNA.
- ✓ 3. Add 1/10 volume 3 M sodium acetate, pH 7.0 and 2.5 volumes cold 100% ethanol.
- x Vortex and incubate <sup>30</sup>~~10~~ min at -80°C.
- ✓ 4. Centrifuge <sup>30</sup>~~10~~ min at 13,000 rpm, 4°C.
- ✓ 5. Discard supernatant (hot!) and add <sup>750 µl</sup> cold 70 % ethanol. Centrifuge <sup>10</sup>~~10~~ min at 13,000 rpm, 4°C. **CHECK!**
- ✓ 6. Discard supernatant (hot!) and dry pellet for 5-10 minutes in the speedvac at 45°C.

+32P  
049274



7. Resuspend pellet in 30  $\mu$ l hybridisation buffer (20 mM HEPES, 0.4 M NaCl, 80% formamide). Vortex for <sup>15 sec</sup> 5 min, then incubate 5 min at 50°C, vortex and spin briefly.
8. Incubate 15 min at 75°C, then 3 h at 50°C to anneal primer (see Note 2).
9. Add 75  $\mu$ l cold 100 % ethanol, vortex and incubate overnight at -80°C. -20 min

## Day 2

1. Precool centrifuge to 4°C and set heating blocks to 37°C and 72°C.
2. Pour sequencing gel. 20
3. Centrifuge annealed primer/RNA sample (stored overnight at -80°C) for 10 min at 13,000 rpm, 4°C.
4. Remove supernatant and add <sup>750  $\mu$ l</sup> 1 ml cold 70 % ethanol. Centrifuge 5 min at 13,000 rpm, 4°C. 10
5. Remove supernatant and dry pellet (approx 20 min).
6. Resuspend pellet in 31  $\mu$ l DEPC-treated water, then add:

10  $\mu$ l 5 x reverse transcriptase buffer (Promega)

1  $\mu$ l 50 mM DTT

5  $\mu$ l 10 mM dNTPs

2.5  $\mu$ l AMV reverse transcriptase (Promega)(see Note 3)

0.6  $\mu$ l RNasin (Promega)

10 $\mu$ l	Go script
5 $\mu$ l	10 mM dNTPs
1 $\mu$ l	Go script
0.6 $\mu$ l	RNasin
MgCl <sub>2</sub>	2 $\mu$ l

7. Incubate 1 h at 37°C. Meanwhile, pre-run sequencing gel.
8. Incubate 10 min at 72°C to inactivate enzyme, then spin briefly.
9. Add 1  $\mu$ l 10 mg/ml RNase and incubate 30 min at 37°C.
10. Add <sup>6.7 3</sup> 6.7  $\mu$ l 4M ~~ammonium~~ <sup>Sodium</sup> acetate, pH 4.8 and 125  $\mu$ l cold 100% ethanol. Centrifuge 10 min, 13,000 rpm, 4°C. 10
11. Remove supernatant and add 1 ml cold 70% ethanol. Centrifuge 7 min at 13,000 rpm, 4°C.
12. Remove supernatant and dry pellet (approx 20 min).
13. Resuspend pellet in 4  $\mu$ l stop solution from T7 sequencing kit.
14. Load 3  $\mu$ l each primer extension reaction on a 6% polyacrylamide sequence gel, along with sequence reactions. Run gel at 60 Watts for approx 2 hours (until dark blue runs off bottom of gel, then a further 3-4 cm).

INCUBATE  
-80°C 30 min.  
-overnight

⇒ Fix gel for 15 min