

In vitro transcriptions assays

Protocol for multiround transcription assays

The final volume of the reaction is 25 μ l – 24 μ l of cocktail and the reaction is started by the addition of 1 μ l polymerase.

- i) Cocktail: for 15 x 24 μ l where RNAP (1 μ l of each in 1 x transcription buffer) is added individually to each reaction:

36 μ l	10 x transcription buffer (Tris/acetate buffer)	2.4y
3.6 μ l	10 mg/ml BSA	0.24y
15 μ l	20 ng/ μ l pSR/lacUV5(short)	y
15 μ l	5 mM GTP	y
15 μ l	5 mM ATP	y
15 μ l	5 mM CTP	y
15 μ l	0.25 mM UTP	y
3.75 μ l	[α^{32} P]-UTP (2.5 μ Ci per reaction)	0.25y
<u>241.65 μl</u>	sterile, RNase-free water	
360 μ l		

- iv) Prepare RNAP dilutions - dilute RNAPs in 1 x transcription buffer + 100 μ g/ml BSA:

Usually use three different concentrations for each polymerase – 4ng/ml, 2ng/ml + 1ng/ml

Reaction number	Cocktail	RNAP dilution
1	24 μ l	1 μ l 1ng/ml
2	24 μ l	1 μ l 2ng/ml
3	24 μ l	1 μ l 4ng/ml
etc		

- iv) Add 1 μ l of RNAP solution to each reaction, each one staggered by 30 seconds, mix by pipetting up and down and moving the tip around, incubate for 15 minutes at 30°C.
- v) Add 25 μ l of stop solution - mix by gentle pipetting.
- vi) Load 15 μ l of each reaction onto the denaturing gel.
- vii) Run gel for ~ 1 hour at 224 V (~ 15V/cm) then dry without fixing.
- viii) Cover gel with cling film and exposed to phosphorimager.
- ix) Analyse the results, measuring the amount of RNA I and the lacUV5 transcripts produced in each lane.

Notes

Transcription buffer:

40 mM Tris acetate pH 8
100 mM KCl
10 mM MgCl₂
1 mM DTT

DNA Template:

PSR is approx. 4000bp \therefore 20ng per reaction is a final conc. Of approx. 0.3 nM.

Protein Concentration:

The DNA concentration is approx. 0.3nM and there are more than 1 promoter per plasmid. If using the polymerase concentrations of 4 μ g/ml, 2 μ g/ml and 1 μ g/ml these equate to 89 μ M, 44 μ M and 22 μ M. This means that the final concentrations in the reactions are 3.52nM, 1.76nM and 0.88nM respectively.